



TracLiNXX



Templates help ensure proper horizontal alignment of tracks during installation

The new **TracLiNXX Track** acts as a spacer, filling the gap between two regular width iNTRAXX tracks and linking the upper and lower tracks together. TracLiNXX is compatible with all iNTRAXX tracks but usability is best when paired with the Standard Track.

TracLiNXX makes it easy to create an iNTRAXX wall using just two Standard tracks—one above and below the TracLiNXX. This Standard Track / TracLiNXX / Standard Track configuration can be repeated as needed to expand the wall without limitation.

A wall of iNTRAXX tracks offers greater flexibility to fine-tune equipment placement exactly where you need it.


Features

- Two recessed channels running the length of TracLiNXX can accommodate LED light strips, if desired
- Easy to clean and decontaminate
- Is compatible with all iNTRAXX track profiles
- Best performance is achieved when TracLiNXX is paired with the Standard Track (0480001)
- A Template Kit (0480293) is recommended to assist with track installation. It includes two templates that help ensure proper horizontal alignment of tracks, allowing attached equipment to move smoothly along the iNTRAXX wall.

Product Part #	0480235
Product Description:	TracLiNXX - 1016mm
Overall Size: mm / inches	(L) 1016 x (W) 87.5 x (D) 12.7 mm (L) 40 x (W) 3.4 x (D) 0.5 in.
Weight	1.5 kg / 3.3 lbs
Load Capacity	N/A
Rating	10g, 22g, 26g



Template Kit (0480293) comprises two templates



Glossary

W: Width

D: Depth

L: Length