


FERNO 309-7 BASKET STRETCHER

DHS-309-7



The Ferno 309-7 Basket Stretcher is a 'Stokes' type litter / basket with a high-density polyethylene patient shell supported and surrounded by a powder coated, aluminium frame. The basket stretcher is an economical addition to any industrial, mining or construction rescue cache.

Features

- Monocoque design
- Full length runners beneath the patient shell provide additional support and assists when dragging over uneven terrain
- Fully exposed outer rail provides good handling and attachment options
- Suitable for vertical or horizontal lifting
- Supplied with three two-piece patient restraints with slot through buckles and Ferno Rapid Release tabs
- Made in Australia

Specifications

•	
Manufacturer's code	DHS-309-7
Length	215 cm
Width	63 cm
Height	24 cm
Weight	13.4 kg
Load Rating	300 kg





Accessories

Ferno Adjustable Lifting Bridle	
BRB-TLB	With steel triple lock karabiners
BRB-TLB-TL	With alloy triple lock karabiners
Basket Stretcher Carry Case	
MCC-387A	Carry case for one-piece basket stretcher
Internal Patient Harness	
BRH-BPH	Basket Patient Harness
BRH-BPH-KIT	Basket Patient Harness and Foot Strap
	coming soon - approx. January 2025
Ferno EARS	
FWEEARS	Anti-Rotation Sail
	minimises basket rotation from helicopter downdraft



FERNO BASKET STRETCHER CARRY CASE MCC-387A (1-PIECE) | MCC-387B (SPLIT)

